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## EN Display box for rehabilitation and anti-stress balls

The display box consists of 24 units of rehabilitation and anti-stress balls of different sizes and densities.

### Characteristics:

1 - Application of the anti-stress ball: the hand rehabilitation ball strengthens finger strength, helps prevent carpal tunnel syndrome, helps relieve stress, helps improve blood circulation, etc.

2 - Three levels of density: The display box includes 3 types of balls with different sizes and density to exert different levels of force, depending on the user's need. The level of resistance is reflected on the ball:

\* Slightly soft level - Yellow

\*\* Moderate level - Red

\*\*\* Hard level - Green

3 - High-quality material: Rehabilitation and anti-stress balls made from PU

### Instructions:

- 1- Grasp the ball and exert force towards the wrist.
- 2- Hold the ball with your fingers on a smooth surface and exert force by pressing the ball on the surface.
- 3- Pinch the ball with your fingers as if it were a pincer.
- 4- Extend the thumb over the ball and apply force to the wrist with the fingers.
- 5- Use force with alternating fingers
- 6- Applying force with the palm of the hand on a smooth surface

### Ball Grip(wrist down)

Grip ball in palm. Turn Hand in so it faces down. Squeeze the ball in the palm. Hold and relax, repeat.

### Opposition

Place ball in the palm of the hand. Keep between thumb and finger being exercised. Squeeze thumb and finger together. Hold and relax hand. Repeat.

### Ball Grip (Wrist Up)

Grip ball in palm. Turn hand but so it faces up. Squeeze the ball with your palm. Hold and relax. Repeat.

### Extend Out

Place ball on the table. Put tips on fingers on the ball. Roll the ball outwards on the table. Repeat.

### Pinch

Place ball between the thumb and index finger. Squeeze together. Hold and relax. Repeat.

### Side Squeezer

Place ball between any two fingers together. Hold and relax. Repeat.

### Thumb Extend

Place ball between bent thumb on two fingers of the same hand. Rolling the ball, extend and straighten the thumb. Repeat.

### Finger Bend

Place ball in the palm of the hand with fingers pressed into the ball. Push fingers into ball as if you are bending your fingers. Hold and then relax. Repeat.

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